

## Product datasheet for **BA1090**

### RNA5-8SN2 Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Human High Density Lipoprotein-3 (HDL3) Subfraction
Species:	Human
Expression Host:	Plasma
Concentration:	1.0 mg/ml (BCA with SDS)
Purity:	Sequential Density Ultracentrifugation at density: 1.120-1.210. Devoid of VLDL, IDL, LDL, and HDL2.
Buffer:	100mM Sodium Chloride, 10 mM Tris, 1mM EDTA, pH 7.4, containing 1mM Sodium Azide as preservative.
Preparation:	Liquid purified protein



[View online »](#)