

## Product datasheet for BA219

### Alpha-2-HS-glycoprotein / AHSB Bovine Protein

#### Product data:

|                        |   |
|------------------------|---|
| Product Type:          | Native Proteins   |
| Description:           | Alpha-2-HS-glycoprotein / AHSB bovine protein, 1 g  |
| Species:               | Bovine  |
| Protein Source:        | Serum   |
| Purity:                | >90% Total impurities: ≤ 0.3% Free N-acetyl-neuraminic acid   |
| Buffer:                | State: Lyophilized, purified<br>Preservative: None<br>Stabilizer: None  |
| Reconstitution Method: | Use sterile distilled water to give a solution of 1 mg/ml.  |
| Preparation:           | Lyophilized, purified   |
| Applications:          | ELISA.  |
| Protein Description:   | Bovine Fetuin from fetal calf serum   |
| Storage:               | Store lyophilized at 2-8°C for 6 months or at -20°C long term.<br>After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term.<br>Avoid repeated freezing and thawing. |
| Stability:             | Shelf life: one year from despatch.   |
| RefSeq:                | <a href="#">NP_776409</a>   |
| Locus ID:              | 280988  |
| Summary:               | alpha-2-HS-glycoprotein   |



[View online »](#)