

Product datasheet for BA350

Mouse IgG3 Mouse Protein

Product data:

| | |
|----------------------|--|
| Product Type: | Native Proteins |
| Description: | Mouse IgG3 mouse protein, 1 mg |
| Species: | Mouse |
| Protein Source: | Cell culture |
| Concentration: | lot specific |
| Purity: | >95% SDS-PAGE >95% |
| Buffer: | State: Liquid Buffer System: Phosphate buffered saline, pH 7.4, 0.09 % Sodium Azide (NaN ₃) |
| Preparation: | Liquid |
| Applications: | ELISA. Control for most immunological assays. |
| Protein Description: | Mouse IgG3 with a kappa light chain |
| Storage: | Store the antibody at 2 - 8 °C up to one month or in aliquots at -20 °C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)