

Product datasheet for BA360

Mouse IgM Mouse Protein

Product data:

Product Type:	Native Proteins
Description:	Mouse IgM mouse protein, 1 mg
Species:	Mouse
Protein Source:	Cell culture
Concentration:	lot specific
Purity:	>95% IEP >95%
Buffer:	State: Liquid Ascites Buffer System: PBS, pH 7.2 containing 0.09% Sodium Azide
Preparation:	Liquid Ascites
Applications:	Western Blot.
Protein Description:	Mouse Immunoglobulin M (IgM) from the MOPC-104E tumour line.
Storage:	Store the protein undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)