

Product datasheet for BA683

RNA5-8SN2 Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Purified Human IgG4. Specificity: A single arc was observed in Immunoelectrophoresis against antisera to whole human serum and human lambda light chain.
Species:	Human
Expression Host:	Myeloma
Concentration:	lot specific
Purity:	> 95 % pure by SDS-PAGE. Preparation: Affinity Chromatography.
Buffer:	PBS containing 0.09% Sodium Azide as preservative.
Preparation:	Liquid purified immunoglobulin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)