

## Product datasheet for **BIN047**

### Hepatitis B pre-S1 Antigen Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Hepatitis B pre-S1 Antigen protein, 0.1 mg
Expression Host:	E. coli
Concentration:	lot specific
Purity:	~95% pure (SDS-PAGE)
Buffer:	Presentation State: Purified State: Liquid Buffer System: PBS, pH 7.2, no preservative
Preparation:	Liquid
Applications:	Suitable for use in ELISA and Western blot immunoassay.
Protein Description:	Pre-S1 antigen, 108 a.a., no Pre-S2. Does not contain a fusion partner.
Storage:	Aliquot and store at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: Six months from despatch.
Synonyms:	HBV preS1
Protein Families:	Suitable for use in ELISA and Western blot immunoassay.



[View online »](#)