

## Product datasheet for **BIN116**

### Coxsackie Virus A9 (P.B. Bozek Strain) Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Coxsackie Virus A9 (P.B. Bozek Strain) protein, 1 ml
Expression Host:	A549
Purity:	Ultra-centrifugation to remove cellular components with the resulting supernatant supplied as antigen.
Buffer:	Presentation State: Lysate State: Liquid lysate Buffer System: No buffer or preservatives.
Preparation:	Liquid lysate
Protein Description:	Coxsackie A9 Antigen, Strain P.B. Bozek
Note:	Internal Thread Vial. Centrifuge Product Before Opening!
	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Store the antigen at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid multiple freeze/thaw cycles.
Stability:	Shelf life: six months from despatch.



[View online »](#)