

Product datasheet for **BIN184**

Influenza A H3N2 (Strain A/Panama/2007/99) Protein

Product data:

Product Type:	Native Proteins
Description:	Influenza A H3N2 Strain A/Panama/2007/99 protein, 0.5 mg
Protein Source:	Chicken eggs
Concentration:	lot specific
Purity:	>90% Ultracentrifugation using 10-40% sucrose gradient.
Buffer:	Presentation State: Purified State: Liquid purified fraction (>90% pure by 7.5% SDS-PAGE). Buffer System: 0.05 M Tris-HCl, pH 8.0 containing 0.1 M Sodium Chloride, 5 mM EDTA with 0.09% Sodium Azide and 0.005% Merthiolate as preservatives.
Preparation:	Liquid purified fraction (>90% pure by 7.5% SDS-PAGE).
Applications:	Serological studies of influenza A virus, immunogen for antibody production. Has been tested with MAb to Influenza A in ELISA.
Protein Description:	Influenza A (A/Panama/2007/99) (H3N2) Influenza Virus Type A (H3N2), Strain A/Panama/2007/99
Note:	Caution: This product has been treated in a manner consistent with methods of inactivation. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Store the antigen at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid multiple freeze/thaw cycles.
Stability:	Shelf life: six months from despatch.
Synonyms:	Seasonal Flu H3N2
Protein Families:	Serological studies of influenza A virus, immunogen for antibody production. Has been tested with MAb to Influenza A in ELISA.



[View online »](#)