

Product datasheet for **BM1093**

E. coli K88a Mouse Monoclonal Antibody [Clone ID: 1050/66]

Product data:

Product Type:	Primary Antibodies
Clone Name:	1050/66
Applications:	ELISA, IF
Recommended Dilution:	ELISA. Immunofluorescence.
Reactivity:	Escherichia coli
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	NCTC 10964 (Nobivac)
Specificity:	In a simple ELISA, using Biogenesis test groups, this antibody is only reactive with <i>E. coli</i> K88a.
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Escherichia coli, EHEC, EPEC, ETEC



[View online »](#)