

Product datasheet for **BM1115**

Feline Immunodeficiency Virus (FIV) gp41 Mouse Monoclonal Antibody [Clone ID: 5A1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	5A1
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western blot. Gives bands at 140 and 41kD in reduced conditions
Reactivity:	Feline Immunodeficiency Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native, FIV extracts.
Specificity:	This antibody reacts with feline immunodeficiency virus / FIV transmembrane protein gp41. It does not cross react with feline leukemia virus, feline herpes virus type 1, feline coronavirus or feline calicivirus.
Formulation:	PBS, containing 0.09% Sodium azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)