

## Product datasheet for **BM1130**

### Hepatitis B Core Antigen / HBcAg Mouse Monoclonal Antibody [Clone ID: H50]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	H50
Applications:	IF
Recommended Dilution:	<b>Immunofluorescence.</b>
Reactivity:	Hepatitis B Virus
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native.
Specificity:	Recognizes the core antigen.
Formulation:	PBS buffer, pH 7.2, 0.05% Sodium Azide as preservative State: Ig Fraction State: Liquid purified IgG bfraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HBV Capsid protein, HBV Core protein, p21.5



[View online »](#)