

## Product datasheet for **BM1160**

### Herpes Virus Type 6 / HHV6 (B-Variant specific) Mouse Monoclonal Antibody [Clone ID: C3108-103]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	C3108-103
Applications:	IF, IHC
Recommended Dilution:	<b>Immunofluorescence:</b> 1/25-1/100. <b>Immunohistochemistry on Paraffin Sections.</b>
Reactivity:	Human Herpesvirus 6
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Purified nucleocapsids.
Specificity:	Clone C3108-103 is specific for the Human Herpes Virus 6, B variant, recognising a viron protein of 101 kDa. Gives a cytoplasmic speckled pattern in IF (See <i>Yamamoto et al.</i> for details). This antibody does not cross react with Human Herpes Virus 6, A variant.
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid purified IgG fraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	HHV-6



[View online »](#)