

## Product datasheet for **BM1251**

### **Pseudomonas aeruginosa serotype 9 Mouse Monoclonal Antibody [Clone ID: 95/159]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	95/159
<b>Applications:</b>	ELISA, IF
<b>Recommended Dilution:</b>	ELISA. Immunofluorescence.
<b>Reactivity:</b>	<i>Pseudomonas aeruginosa</i>
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	<i>P. aeruginosa</i> , serotype 9 (NCTC 11448, ATCC 33356)
<b>Specificity:</b>	In a simple ELISA using Biogenesis Test Groups this antibody is reactive with <i>P. aeruginosa</i> , serotype 9.
<b>Formulation:</b>	PBS buffer, pH 7.2, 0.1% Sodium Azide State: Purified State: Liquid purified IgG
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Protein A chromatography
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	The antibody can be shipped at 2-8°C. Store (in aliquots) at -20°C only. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)