

## Product datasheet for **BM1302**

### Schistosoma mansoni Mouse Monoclonal Antibody [Clone ID: BGN/1A7]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	BGN/1A7
Applications:	ELISA, IF
Recommended Dilution:	<b>ELISA.</b> <b>Immunofluorescence.</b>
Reactivity:	Schistosoma mansoni
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Native <i>S. mansoni cercariae</i> .
Specificity:	Reacts with antigen present in the urine of patients with Schistosomiasis. Reacts with a TCA soluble adult worm antigen in IF/ELISA.
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	<i>S. mansoni</i>



[View online »](#)