

## Product datasheet for **BM1315**

### Shigella flexneri (Serotype 6) Mouse Monoclonal Antibody [Clone ID: 306/305]

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Clone Name:           | 306/305  |
| Applications:         | ELISA, IF  |
| Recommended Dilution: | ELISA.   |
| Reactivity:           | Vibrio cholerae  |
| Host:                 | Mouse  |
| Isotype:              | IgG  |
| Clonality:            | Monoclonal   |
| Immunogen:            | <i>S. flexneri</i> , NCTC 9729   |
| Specificity:          | This antibody is specific for <i>Shigella flexneri</i> serotype 6. <i>S. flexneri</i> is a gram negative, rod shaped bacteria that is very closely related to <i>E. coli</i> and <i>Salmonella</i> . |
| Formulation:          | PBS, pH 7.4 containing 0.09% Sodium Azide<br>State: Purified<br>State: Liquid purified IgG   |
| Concentration:        | lot specific   |
| Purification:         | Affinity chromatography on Protein G   |
| Conjugation:          | Unconjugated   |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing.  |
| Stability:            | Shelf life: one year from despatch.  |



[View online »](#)