

Product datasheet for **BM1320**

Staphylococcus aureus / MRSA Mouse Monoclonal Antibody [Clone ID: 332/423]

Product data:

Product Type:	Primary Antibodies
Clone Name:	332/423
Applications:	ELISA, IF, WB
Recommended Dilution:	ELISA. Western Blot. Immunofluorescence. Reactive in immunohistology with isolates fixed in Phenol, Formalin and by boiling.
Reactivity:	Staphylococcus aureus
Host:	Mouse
Isotype:	IgG2b
Clonality:	Monoclonal
Immunogen:	MRSA
Specificity:	Recognises methicillin-resistant Staphylococcus aureus (MRSA), also known as oxacillin-resistant Staphylococcus aureus (ORSA).
Formulation:	PBS, pH 7.4, containing 0.09% sodium azide as preservative State: Purified State: Liquid IgG fraction
Concentration:	lot specific
Purification:	Protein G affinity chromatography
Conjugation:	Unconjugated
Storage:	Store at 2-8°C or at -20°C if preferred. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use.
Stability:	Shelf life: One year from despatch.



[View online »](#)