

Product datasheet for **BM1350**

VC1456 Mouse Monoclonal Antibody [Clone ID: 2/63]

Product data:

Product Type:	Primary Antibodies
Clone Name:	2/63
Applications:	ELISA, IHC
Recommended Dilution:	ELISA, Immunohistochemistry (Formaldehyde is recommended as fixative).
Reactivity:	Vibrio cholerae
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Purified native cholera Toxin beta
Specificity:	Reacts specifically with beta chain of Vibrio cholerae toxin in a simple ELISA.
Formulation:	PBS, pH 7.4, containing 0.09% sodium azide. State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store at 2-8°C for up to one month or aliquot and store at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P01556
Synonyms:	Cholera enterotoxin B chain, Cholera enterotoxin gamma chain, Cholera toxin B subunit, CTxB, CTxB, VC_1456



[View online »](#)