

Product datasheet for **BM2107**

Benzoyllecgonine Mouse Monoclonal Antibody [Clone ID: IP3G2]

Product data:

Product Type:	Primary Antibodies
Clone Name:	IP3G2
Applications:	ELISA, R
Recommended Dilution:	For competitive RIA or EIA for BE or Cocaine quantitation.
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Benzoyllecgonine - KLH
Specificity:	This antibody detects Benzoyllecgonine. <i>Cross reactivities:</i> Cocaine: > 100% Ecgonine: 3.2% Ecgonine Methylester: < 0.1%
Formulation:	PBS, pH 7.2 State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	DEAE Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Benzoyllecgonine is the predominant cocaine metabolite. This compound is not biologically active but is used as an indicator of cocaine use.
Synonyms:	Cocaine



[View online »](#)