

Product data:

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for BM2226P

IGF1 Receptor (IGF1R) (alpha chain) Mouse Monoclonal Antibody [Clone ID: 1H7]

Product Type: Primary Antibodies Clone Name: 1H7 FC, IHC, WB **Applications:** Recommended Dilution: Western Blot. Immunohistochemistry on Paraffin Sections: 1/10-1/40. Requires Protein digestion pre-treatment of paraffin sections (e.g. trypsin or pronase). **Flow Cytometry**: Use 10 µl of 1/10 diluted antibody to label 16 cells in 100 µl. Recommended Positive Control: Pancreas, placenta. **Reactivity:** Human Host: Mouse Isotype: lgG1 **Clonality:** Monoclonal Immunogen: Purified Human placental IGF-I receptor. Specificity: This antibody recognizes CD221, a 155kD receptor tyrosine kinase, also known as Insulin-like growth factor I receptor (IGF-I Receptor). CD221 is composed of two extracellular alphasubunits and two transmembrane betasubunits. Clone 1H7 recognizes an epitope in the alpha subunits of CD221, demonstrated by Western blotting. CD221 is expressed in a range of tissues, including kidney, liver, placenta, mammary gland, brain, ovary and skin. The ligands for CD221 include IGF-I and IGF-II, which bind to CD221 and activate tyrosine kinase activity, resulting in phosphorylation of several intracellular signalling proteins. Clone 1H7 is reported to partically block binding of IGF-I and IGF-II to CD221. We recommend the use of BM2226LE for this purpose.

Formulation:PBS
State: Purified
State: Liquid purified IgG fraction from tissue culture supernatant
Preservative: 0.09% Sodium AzideConcentration:lot specificPurification:Affinity Chromatography on Protein G



/iew online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US


	IGF1 Receptor (IGF1R) (alpha chain) Mouse Monoclonal Antibody [Clone ID: 1H7] – BM2226P
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	insulin like growth factor 1 receptor
Database Link:	Entrez Gene 3480 Human P08069
Background:	There are two kinds of IGF receptors based on their relative affinities for IGF1 versus IGF2. The IGF1 receptor (also known as CD221) prefers IGF1 over IGF2 and weakly binds insulin. IGF1 receptor is a disulfide-linked heterotetrameric transmembrane protein consisting of two alpha (130 kD) and two beta (95 kD) subunits. Both the alpha and beta subunits are encoded within a single receptor precursor cDNA. The IGF1 receptor is therefore similar in structure to the insulin receptor. The proreceptor polypeptide is proteolytically cleaved and disulfide- linked to yield the mature heterotetrameric receptor. The IGF1 receptor is highly expressed in all cell types and tissues and is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival.
Synonyms:	IGF-I receptor, IGF1 Receptor

Product images:

555 416-502277-138-0-10¹⁰ 10¹ 10² 10³ 10⁴

Staining of human peripheral blood granulocytes with mouse anti human CD221

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Staining of human peripheral blood ganulocytes with mouse anti human CD221: Alexa FluorR 488

Staining of human peripheral blood ganulocytes with mouse anti human CD221: Low Endotoxin

Staining of human peripheral blood granulocytes with mouse anti human CD211:RPE

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US