

Product datasheet for **BM2431P**

CD44 Mouse Monoclonal Antibody [Clone ID: 156-3C11]


Product data:

Product Type:	Primary Antibodies
Clone Name:	156-3C11
Applications:	FC, IHC, IP, WB
Recommended Dilution:	Flow Cytometry: Use 10 µl of 1/50-1/100 diluted CD44 antibody to label 1x10 ⁶ cells in 100 µl. Immunoprecipitation. Western Blot: Detects a band of approximately 80kDa in HeLa cell lysates. Immunohistochemistry on Frozen and Paraffin Sections: 1/100. This product requires antigen retrieval using heat treatment prior to staining of paraffin sections. Sodium citrate buffer pH 6.0 is recommended for this purpose. Recommended Positive Control: Human tonsil.
Reactivity:	Human, Monkey
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Stimulated Human leucocytes.
Specificity:	This antibody recognises CD44, a 90kD type I transmembrane protein which is otherwise known as HCAM or CD44s. Various isoforms of CD44 exist. The 90kD standard form, recognised by this antibody, is expressed on leucocytes, erythrocytes, white matter of the brain and some epithelial cells of the breast and small intestine. Antibodies produced by clone 156-3C11 recognise epitope 3, which is a protease resistant epitope on the CD44 molecule.
Formulation:	PBS, pH 7.4 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein G
Conjugation:	Unconjugated



[View online »](#)

Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD44 molecule (Indian blood group)
Database Link:	Entrez Gene 960 Human P16070
Background:	CD44 is a receptor for hyaluronic acid (HA) and is involved in cell-cell interactions, cell adhesion and migration. CD44 also participates in a wide variety of cellular functions including lymphocyte activation, recirculation and homing. CD44 expression may be up-regulated upon some carcinomas, and it has been speculated that this may be related to metastatic potential.
Synonyms:	LHR, MDU2, MDU3, MIC4, CDw44, Epican, ECMR-III, HUTCH-I, Heparan sulfate proteoglycan, Hermes antigen, Hyaluronate receptor, PGP-1

Product images:

Western blot analysis of HeLa whole cell lysate probed with Mouse anti Human CD44 antibody (M2431P) followed by HRP conjugated Goat anti Mouse IgG, visualized by chemiluminescence.