

Product datasheet for **BM2711**

Gentamycin Mouse Monoclonal Antibody [Clone ID: 102]

Product data:

Product Type:	Primary Antibodies
Clone Name:	102
Applications:	ELISA
Recommended Dilution:	May be used in Immunoassays.
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	Gentamicin conjugated to KLH
Specificity:	This antibody specifically recognizes Gentamicin.
Formulation:	PBS, pH 7.4 without preservatives. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Purification:	Protein G chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Gentamycin is an aminoglycoside antibiotic used to treat infections caused by many different types of bacteria. Gentamycin is usually administered by intravenous infusion or intramuscular injection.



[View online »](#)