

## Product datasheet for **BM3039**

### **Yersinia enterocolitica 0:9 Mouse Monoclonal Antibody [Clone ID: 8E9]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	8E9
<b>Applications:</b>	AGG, IF
<b>Recommended Dilution:</b>	Immunofluorescence: 1:20. Agglutination: Ready-to-use.
<b>Reactivity:</b>	Yersinia enterocolitica
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG3
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Yersinia enterocolitica serogroup O:9 strain
<b>Specificity:</b>	8E9 specifically reacts with Yersinia enterocolitica serogroup O:9 strains in both agglutination and immunofluorescence assays. The antibody does not react with Y. enterocolitica O:3 and O:8 serogroups.
<b>Formulation:</b>	PBS with 0.5% BSA and 0.09% Sodium Azide State: Purified State: Liquid purified IgG fraction
<b>Purification:</b>	Protein-A affinity-purified
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)