

## Product datasheet for **BM3055**

### **E. coli (mixed serotypes) Mouse Monoclonal Antibody [Clone ID: 405]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	405
<b>Applications:</b>	ELISA, ID
<b>Recommended Dilution:</b>	Indirect ELISA. RID.
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Pool of four E. coli serotypes, which are common to human E. coli diseases. These serotypes are 018:B21, 044:K74, 0112:B11, 0125:B15.
<b>Specificity:</b>	Specific for E. coli serotypes O2a/2b, O7, O11, O18, O44, O75, O122, and O125. Does not cross-react with E. aerogenes, P. vulgaris, K. pneumoniae, P. mirabilis, P. vulgaris, Ps. Aeruginosa or S. marcescens.
<b>Formulation:</b>	State: Ascites State: Liquid ascites with 0.09% sodium azide as preservative.
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody (in aliquots) at -20°C . Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Synonyms:</b>	Escherichia coli, EHEC, EPEC, ETEC



[View online »](#)