

Product datasheet for **BM3074P**

Streptococcus Group A Mouse Monoclonal Antibody [Clone ID: 501]

Product data:

Product Type:	Primary Antibodies
Clone Name:	501
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Streptococcus sp.
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Immunogen:	UV-inactivated Group A Streptococcus cells
Specificity:	This antibody reacts to Streptococcus Group A. Does not cross react with Streptococcus groups B-G, <i>Streptococcus pneumoniae</i> , <i>Staphylococcus aureus</i> , <i>Corynebacterium diphtheriae</i> , <i>Streptococcus salivarius</i> or <i>Micrococcus mucilagines</i> .
Formulation:	PBS, pH 7.4 containing 0.09% Sodium azide as preservative State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C to -70°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)