

## Product datasheet for **BM3204**

### Measles virus Mouse Monoclonal Antibody [Clone ID: 166]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	166
Applications:	ELISA, IF
Recommended Dilution:	<b>ELISA</b> <b>Immunofluorescence.</b>
Host:	Mouse
Isotype:	IgG
Clonality:	Monoclonal
Immunogen:	Whole virus, Edmunston strain. Source: Ascites.
Specificity:	This antibody is specific for Measles HA.
Formulation:	0.01 M PBS, pH 7.2 State: Purified State: Liquid purified Ig fraction (>90% pure) Stabilizer: None Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Rubeola virus



[View online »](#)