

Product datasheet for **BM3353**

Yersinia pestis V Antigen Mouse Monoclonal Antibody [Clone ID: Va48]

Product data:

Product Type:	Primary Antibodies
Clone Name:	Va48
Applications:	ELISA, WB
Recommended Dilution:	ELISA: <u>Recommended Pairs for Sandwich Immunoassay:</u> Detection/Capture: BM3353/BM3351 Western blot.
Reactivity:	Yersinia pestis
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Recombinant full-length V antigen of Yersinia pestis
Specificity:	Yersinia pestis V antigen.
Formulation:	PBS, pH 7.4 with 0.09% Sodium Azide as preservative. State: Purified State: Liquid purified Ig fraction (>95% pure).
Concentration:	lot specific
Purification:	Protein G Chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C
Stability:	Shelf life: one year from despatch.



[View online »](#)