

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for BM4041F

CD163 Mouse Monoclonal Antibody [Clone ID: 5C6-FAT]

Product data:

Product Type:	Primary Antibodies
Clone Name:	5C6-FAT
Applications:	FC
Recommended Dilution:	Flow Cytometry: Use 50 μl of 1/10-1/20 diluted antibody to label 5x1 ^g cells. <i>Suggested Positive Control:</i> Human monocytes.
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Human CD163.
Specificity:	 Monoclonal antibody 5C6 FAT recognizes a membrane Glycoprotein on Monocytes and Macrophages which is expressed in intermediate and late inflammatory stages. Antigen Distribution Isolated Cells: Monocytes, particularly after dexamethasone treatment or after 2-5 days in culture. Does not react with lymphocytes, granulocytes or platelets. Tissue Sections: Positive staining can be observed in the skin (histiocytes), gut, Kupffer cells, few alveolar macrophages, a major population of macrophages in the placenta, varying degrees of macrophages in inflamed tissues, including tumorous tissue depending on the inflammatory stage. Red pulp, but not white pulp macrophages of the spleen, and cortical macrophages of the thymus are detected. Macrophages in the synovialis of patients with rheumatoid arthritis. In alveolar macrophages and in Kupffer cells a double staining can be observed with monoclonal antibody 25F9 (product <i>CatNo</i> BM4022) which is not the case in other tissues. The antigen is also found abundantly in human term placenta, and regularly in acute and chronic inflammatory lesions.



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US CD163 Mouse Monoclonal Antibody [Clone ID: 5C6-FAT] – BM4041F

Formulation:	PBS, pH 7.2 Label: FITC State: Liquid purified IgG fraction Stabilizer: 10 mg/ml BSA Preservative: 0.05% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	0.2 mg/ml
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	CD163 molecule
Database Link:	<u>Entrez Gene 9332 Human</u> <u>Q86VB7</u>
Background:	CD163 is a scavenger receptor for the haemoglobin-haptoglobin complex, and is upregulated by glucocorticoids and IL-10. The extracellular portion of the receptor is regularly shed and can be found in the circulation. An important function of CD163 seems to be in the adhesion of monocytes to activated endothelial cells. CD163-positive cells include skin histiocytes, Kupffer cells, spleen macrophages of the red pulp, and some thymus macrophages. The antigen is also found abundantly in human term placenta, and regularly in acute and chronic inflammatory lesions.
Synonyms:	M130, Hemoglobin scavenger receptor, Macrophage marker

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2024 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US