

### Product datasheet for BM4041X

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# CD163 Mouse Monoclonal Antibody [Clone ID: 5C6-FAT]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: 5C6-FAT

**Applications:** ELISA, FC, IHC, WB

Recommended Dilution: ELISA.

**Immunohistochemistry on Frozen Sections:** 0.5 μg/ml (1/4000)

Immunohistochemistry on Paraffin Sections: 1 µg/ml (1/2000) (Proteinase K pretreatment

for antigen retrieval is recommended). *Suggested Positive Control:* Human placenta.

Has been described to work in FACS and Western blots.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal Human CD163.

The epitope has not been further characterized.

**Specificity:** Monoclonal antibody 5C6 FAT recognizes a membrane glycoprotein on monocytes and

macrophages which is expressed in intermediate and late inflammatory stages.

Antigen distribution:

**Isolated cells:** Monocytes, particularly after dexamethasone treatment or after 2-5 days in

culture. Does not react with lymphocytes, granulocytes or platelets.

**Tissue Sections:** Positive staining can be observed in the skin (histiocytes), gut, Kupffer cells, few alveolar macrophages, a major population of macrophages in the placenta, varying degrees of macrophages in inflamed tissues, including tumorous tissue depending on the inflammatory stage. Red pulp, but not white pulp macrophages of the spleen, and cortical macrophages of the thymus are detected. Macrophages in the synovialis of patients with rheumatoid arthritis. In alveolar macrophages and in Kupffer cells a double staining can be observed with monoclonal antibody 25F9 (product BM4022) which is not the case in other

tissues.





### CD163 Mouse Monoclonal Antibody [Clone ID: 5C6-FAT] - BM4041X

**Formulation:** PBS, pH 7.2

State: Purified

State: Liquid purified IgG fraction.

Concentration: 2.0 mg/ml

**Purification:** Affinity Chromatography.

Conjugation: Unconjugated

**Storage:** Store original vial at -70°C.

Aliquots of stock solution: at -70°C. Do not freeze working dilutions. Avoid repeated freezing and thawing.

**Stability:** Shelf life: one year from despatch.

Gene Name: CD163 molecule

Database Link: Entrez Gene 9332 Human

Q86VB7


**Background:** CD163 is a scavenger receptor for the haemoglobin-haptoglobin complex, and is upregulated

by glucocorticoids and IL-10. The extracellular portion of the receptor is regularly shed and can be found in the circulation. An important function of CD163 seems to be in the adhesion of monocytes to activated endothelial cells. CD163-positive cells include skin histiocytes, Kupffer cells, spleen macrophages of the red pulp, and some thymus macrophages. The antigen is also found abundantly in human term placenta, and regularly in acute and chronic

inflammatory lesions.


Synonyms: M130, Hemoglobin scavenger receptor, Macrophage marker

## **Product images:**




Paraffin Sections of Human Lung using 5C6-FAT






Paraffin Sections of Human Placenta using 5C6-FAT



Frozen Sections of Human Placenta using 5C6-FAT



Paraffin Sections of Human Tonsils using 5C6-FAT