

Product datasheet for **BM469**

Human IgA (Fc specific) Mouse Monoclonal Antibody [Clone ID: 15D6]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA (Fc specific) Mouse Monoclonal Antibody [Clone ID: 15D6]
Clone Name:	15D6
Applications:	ELISA
Recommended Dilution:	ELISA: 1/500-1/1000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified Human IgA.
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Protein G Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)