

Product datasheet for **BM472**

Human IgA1 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 3C10]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgA1 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 3C10]
Clone Name:	3C10
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgA, from colostrum
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)