

Product datasheet for **BM474**

Human IgD (Fc specific) Mouse Monoclonal Antibody [Clone ID: LE2]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Human IgD (Fc specific) Mouse Monoclonal Antibody [Clone ID: LE2] |
| Clone Name: | LE2 |
| Applications: | IHC |
| Recommended Dilution: | Flow cytometry. Immunohistochemistry on frozen tissue sections. |
| Reactivity: | Human |
| Host: | Mouse |
| Immunogen: | Partially purified human IgD (Fc) |
| Isotype: | IgG1 |
| Formulation: | State: Ig Fraction State: Liquid purified IgG containing 0.09% Sodium Azide |
| Purification: | Affinity chromatography |
| Conjugation: | Unconjugated |
| Storage: | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)