

Product datasheet for **BM480**

Human IgG2 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: 3C7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG2 (Hinge Region) Mouse Monoclonal Antibody [Clone ID: 3C7]
Clone Name:	3C7
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgG from pooled serum
Isotype:	IgG1
Formulation:	Phosphate Buffered Saline State: Purified State: Liquid purified Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)