

Product datasheet for **BM489**

Human IgM Mouse Monoclonal Antibody [Clone ID: 47.2.E9]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Mouse Monoclonal Antibody [Clone ID: 47.2.E9]
Clone Name:	47.2.E9
Applications:	ELISA
Recommended Dilution:	ELISA (1/5000-1/10000)
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgM.
Isotype:	IgG
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)