

Product datasheet for **BM490**

Human IgM Rat Monoclonal Antibody [Clone ID: HM-7]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM Rat Monoclonal Antibody [Clone ID: HM-7]
Clone Name:	HM-7
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Rat
Immunogen:	Purified Human IgM.
Isotype:	IgM
Formulation:	PBS State: Aff - Purified State: Liquid purified IgM fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)