

Product datasheet for **BM498**

Rat IgA Mouse Monoclonal Antibody [Clone ID: MARA-2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgA Mouse Monoclonal Antibody [Clone ID: MARA-2]
Clone Name:	MARA-2
Applications:	ELISA, R
Recommended Dilution:	ELISA. Radioimmunoassays.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Native IgA from Rat (<i>Rattus norvegicus</i>).
Isotype:	IgG1
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Stabilizer: 50% Glycerol Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)