

Product datasheet for **BM499**

Rabbit IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: RG-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG (heavy chain) Rat Monoclonal Antibody [Clone ID: RG-1]
Clone Name:	RG-1
Applications:	ELISA, IF, IHC
Recommended Dilution:	ELISA. Immunofluorescence.
Reactivity:	Rabbit
Host:	Rat
Immunogen:	Purified Rabbit immunoglobulin gamma.
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide State: Aff - Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)