

Product datasheet for **BM582**

GNRHR Mouse Monoclonal Antibody [Clone ID: A9E4]

Product data:

Product Type:	Primary Antibodies
Clone Name:	A9E4
Applications:	ELISA, FC, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Flow Cytometry. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Synthetic peptide corresponding to the amino acids 1-29 of the Human GnRH receptor extracellular domain conjugated to BSA.
Specificity:	This Monoclonal antibody <i>A9E4</i> reacts with GnRH receptors in the anterior pituitary.
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	gonadotropin releasing hormone receptor
Database Link:	Entrez Gene 2798 Human P30968



[View online »](#)

Background:

Gonadotropin releasing hormone (GnRH) receptor is a protein of 54-60 kDa. GnRH receptor contains seven hydrophobic transmembrane domains connected by hydrophilic extracellular and intracellular loops characteristic of G-protein coupled receptors. GnRH stimulates the gonadotrophs of the anterior pituitary to secrete luteinizing hormone (LH) as well as follicle-stimulated hormone (FSH). The receptor contains of seven hydrophobic transmembrane domains connected by hydrophilic extracellular, and intracellular loops characteristic of G-protein coupled receptors.

Synonyms:

Gonadotropin-releasing hormone receptor, GnRH-R, GRHR