

Product datasheet for **BM948**

Dengue Virus 2 Mouse Monoclonal Antibody [Clone ID: 3H5.1]

Product data:

Product Type:	Primary Antibodies
Clone Name:	3H5.1
Applications:	IF
Recommended Dilution:	Immunofluorescence: 1/100 - 1/400.
Reactivity:	Dengue Virus Types 2
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	Dengue type 2 antigen (New Guinea C)
Specificity:	This antibody reacts with Dengue virus type 2 recognizes 2645-5057.
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography on Protein G
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Dengue virus type 2, Dengue 2



[View online »](#)