

Product datasheet for **BP2023**

FeLV gag p27 Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Suitable for use in Western blot. Also suitable for conjugation purposes.
Reactivity:	Feline Leukemia Virus
Host:	Goat
Clonality:	Polyclonal
Immunogen:	Purified FeLV p27
Specificity:	FeLV p27 (Western blot). Non-reactive with normal feline serum.
Formulation:	0.01 M PBS, pH 7.2, containing 0.09 % Sodium azide as preservative without stabilizing proteins. State: Ig Fraction State: Liquid purified Ig fraction (>95% pure)
Concentration:	lot specific
Purification:	Sodium sulfate precipitation and ion exchange chromatography.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Feline Leukemia Virus



[View online »](#)