

Product datasheet for **BP406**

FCY1 Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA (1/500 - 1/2,500). Western blot (1/100 - 1/500).
Reactivity:	Yeast
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Highly pure Cytosine Deaminase from yeast
Specificity:	This antibody is specific for cytosine deaminase. It recognises the native antigen in ELISA, and in Western blot two bands were visualised by ECL at 13kD and 17kD.
Formulation:	PBS containing 0.09% sodium azide State: Ig Fraction State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Q12178
Synonyms:	Cytosine aminohydrolase, FCY1



[View online »](#)