

## Product datasheet for BP5024

### Adeno-Associated Virus 2 / AAV2 (VP1 + VP2 + VP3) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IP, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunoprecipitation:</b> 1//10. <b>Western Blot:</b> 1/200-1/500 for Immunoblotting. <b>Immunofluorescence Microscopy:</b> 1/50-1/100. <b>Dilution Buffer:</b> PBS.
Reactivity:	Adeno-Associated Virus 2
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Recombinant adeno-associated virus (AAV-2) capsid proteins, expressed in baculovirus.
Specificity:	The antibody reacts with VP1, VP2 and VP3 of Adenoassociated virus-2.
Formulation:	State: Serum State: Liquid Stabilized antiserum with 0.09% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C.
Stability:	Shelf life: one year from despatch.
Synonyms:	AAV-2

