

## Product datasheet for **BP722**

### Paraquat Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	R
Recommended Dilution:	Radioimmunoassays (has been used in a 125I RIA with derivative of Paraquat).
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Paraquat N-valeric acid derivative (N- (4-carboxybutyl)-N'-methyl-bipyridilium) conjugated to keyhole limpet haemocyanin (KLH).
Specificity:	This antibody is specific for: Paraquat 100 %, Monoquat 5 %, Diquat < 0.1 %, Norfamaquat < 0.1 %.
Formulation:	PBS, pH 7.2, 0.09 % Sodium Azide State: Serum State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Protein G chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20 °C. Can be shipped at 2 - 8 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Synonyms:	Methyl viologen



[View online »](#)