

## Product datasheet for **BP879**

### Vancomycin Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA.
Host:	Sheep
Clonality:	Polyclonal
Immunogen:	Vancomycin conjugated to Keyhole limpet Haemocyanin.
Specificity:	This antibody is specific for Vancomycin.
Formulation:	State: Serum State: Liquid serum without preservatives
Conjugation:	Unconjugated
Storage:	Store (in aliquots) at -20 °C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Vancomycin is a glycopeptide antibiotic that attacks the D-alanyl D-alanine component of the bacterial cell wall. Up until recently, it was used in the prophylaxis and treatment of Gram-positive bacterial infections. Some species have now started to develop resistance due to the structural modification of various surface glycoproteins.



[View online »](#)