

## Product datasheet for **CB500543**

### Frozen Tissue Blocks, Ovary

#### Product data:

|   |                                    |
|---|------------------------------------|
| Product Type:                                     | Frozen Blocks                      |
| Disease State:                                    | Cancer                             |
| Diagnosis Category:                               | Ovarian Cancer                     |
| Tissue:   | Ovary                              |
| Frozen Sample ID:                                 | FR00002C31                         |
| Case ID:  | CI0000000069                       |
| Tissue of Origin:                                 | Ovary                              |
| Site of Finding:                                  | Ovary                              |
| Appearance:                                       | Tumor                              |
| Sample Diagnosis (from Pathology Verification):   | Adenocarcinoma of ovary, serous    |
| Normal:   | 0%                                 |
| Lesional:   | 0%                                 |
| Tumor:  | 65%                                |
| Tumor Hypo/Acellular Stroma:                      | 0%                                 |
| Tumor Hypercellular Stroma:                       | 25%                                |
| Necrosis:   | 10%                                |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of ovary, serous    |
| Age:  | 51                                 |
| Gender:   | Female                             |
| Tumor Grade:                                      | FIGO G2: Moderately differentiated |
| Minimum Stage Grouping:                           | III                                |



[View online »](#)

T (Extent of Primary Tumor): pT3


N (Lymph node metastasis): pNX

M (Distant metastasis): pMX


Storage: Store frozen tissue blocks at -80°C.

Stability: Frozen tissue blocks are stable for at least one year from date of shipping when stored at -80°C.

### Product images:



4x Image



20x Image