

Product datasheet for CB528174

Frozen Tissue Blocks, Tonsil

Product data:

| | |
|---|--------------------------------|
| Product Type: | Frozen Blocks |
| Disease State: | Normal |
| Diagnosis Category: | Normal |
| Tissue: | Tonsil |
| Frozen Sample ID: | FR0002D7CF |
| Case ID: | CI0000010426 |
| Tissue of Origin: | Tonsil |
| Site of Finding: | Tonsil |
| Appearance: | Normal |
| Normal: | 100% |
| Lesional: | 0% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| CASE Diagnosis (from medical center path report): | Lymphoid hyperplasia of tonsil |
| Age: | 34 |
| Gender: | Female |
| Tumor Grade: | n/a |
| Minimum Stage Grouping: | n/a |
| T (Extent of Primary Tumor): | n/a |


[View online »](#)

N (Lymph node metastasis): n/a

M (Distant metastasis): n/a

Storage: Store frozen tissue blocks at -80°C.

Stability: Frozen tissue blocks are stable for at least one year from date of shipping when stored at -80°C.