

Product datasheet for **CB529851**

Frozen Tissue Blocks, Lung

Product data:

| | |
|---|---|
| Product Type: | Frozen Blocks |
| Disease State: | Other Disease |
| Diagnosis Category: | Chronic Obstructive Pulmonary Disease (COPD) |
| Tissue: | Lung |
| Frozen Sample ID: | FR00030FC8 |
| Case ID: | CI0000011577 |
| Tissue of Origin: | Lung |
| Site of Finding: | Lung |
| Appearance: | Lesion |
| Sample Diagnosis (from Pathology Verification): | Emphysema |
| Normal: | 40% |
| Lesional: | 60% |
| Tumor: | 0% |
| Tumor Hypo/Acellular Stroma: | 0% |
| Tumor Hypercellular Stroma: | 0% |
| Necrosis: | 0% |
| Pathology Verification notes from H&E review: | 70% Alveoli, 10% Bronchioles, 20% Fibrovascular septa |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of lung |
| Age: | 67 |
| Gender: | Male |
| Tumor Grade: | n/a |




[View online »](#)

Minimum Stage Grouping: n/a


Storage: Store frozen tissue blocks at -80°C.

Stability: Frozen tissue blocks are stable for at least one year from date of shipping when stored at -80°C.

Product images:



4x Image



20x Image