

## **Product datasheet for CB530799**

## **Frozen Tissue Blocks, Tonsil**

## **Product data:**

**Product Type:** Frozen Blocks **Disease State:** Other Disease **Diagnosis Category:** Hyperplasia

Tissue: Tonsil

Frozen Sample ID: FR00032191 Case ID: CI0000012346

Tissue of Origin: Tonsil Site of Finding: Tonsil Appearance: Lesion

Sample Diagnosis (from Pathology Verification):

Follicular hyperplasia of tonsil

0% Normal: Lesional: 100% 0% Tumor: 0%

Tumor Hypo/Acellular Stroma:

**Tumor Hypercellular** 0%

Stroma:

**Necrosis:** 

**Pathology Verification** notes from H&E review: 5% Squamous epithelium, 95% Lymphoid nodules

**CASE Diagnosis (from** medical center path

report):

Lymphoid hyperplasia of tonsil

44 Age:

Gender: **Female** 

Store frozen tissue blocks at -80°C. Storage:



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn


Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com




Stability:

Frozen tissue blocks are stable for at least one year from date of shipping when stored at  $80^{\circ}$ C.

## **Product images:**



4x Image



20x Image