

## Product datasheet for **CCRT501**

### TissueScan, Crohn's/Colitis cDNA Array I

#### Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/ccrt501.xls
Diagnosis Category:	Crohns and Colitis
Tissue:	Intestine, large and small
Array Name:	Crohn's/Colitis cDNA Array I


#### Product images:



TissueScan product image.




[View online »](#)



A Human Crohn's and Colitis real time PCR Gene Expression panel was used in real-time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed on a PE7700.

Fluorescence reading - cycle number



A Human Crohn's and Colitis real time PCR Gene Expression panel was used in SYBR Green I based real time PCR. The expression profile of human tumor necrosis factor (TNF) gene (NM\_000594) was detected using SYBR Green I protocol and performed on a PE7700.

Normalized relative expression among 48 tissue samples