

Product datasheet for **CD563105**

Tissue Genomic DNA, Liver

Product data:

Product Type:	Genomic DNA
Disease State:	Normal
Diagnosis Category:	Normal
Tissue:	Liver
Frozen Sample ID:	FR5B33B04F
Case ID:	CI0000008106
Prep:	B
Vial ID:	DN00000BD1
Tissue of Origin:	Liver
Site of Finding:	Liver
Appearance:	Normal
Sample Diagnosis (from Pathology Verification):	Within normal limits
Normal:	100%
Lesional:	0%
Tumor:	0%
Tumor Hypo/Acellular Stroma:	0%
Tumor Hypercellular Stroma:	0%
Necrosis:	0%
Pathology Verification notes from H&E review:	Non Tumor Structures: 95% Lobules, 5% Triads
Bioanalyzer Ratio (28S/18S):	2.08


[View online »](#)

CASE Diagnosis (from medical center path report):	Adenocarcinoma of colon, metastatic
Age:	75
Gender:	Male
Race:	White or Caucasian
Tumor Grade:	Not applicable
Minimum Stage Grouping:	Not applicable
T (Extent of Primary Tumor):	Not applicable
N (Lymph node metastasis):	Not applicable
M (Distant metastasis):	Not applicable
Storage:	Store the total DNA -80°C.
Stability:	These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws.