

## Product datasheet for **CD563144**

### Tissue Genomic DNA, Endometrium


#### Product data:

|   |   |
|---|---|
| Product Type:                                     | Genomic DNA                                 |
| Disease State:                                    | Cancer                                      |
| Diagnosis Category:                               | Endometrial Cancer                          |
| Tissue:   | Endometrium                                 |
| Frozen Sample ID:                                 | FR00001047                                  |
| Case ID:  | CU0000000462                                |
| Prep:   | A   |
| Vial ID:  | DN000000BC                                  |
| Tissue of Origin:                                 | Endometrium                                 |
| Site of Finding:                                  | Uterus                                      |
| Appearance:                                       | Tumor                                       |
| Sample Diagnosis (from Pathology Verification):   | Adenocarcinoma of endometrium, endometrioid |
| Normal:   | 0%  |
| Lesional:   | 0%  |
| Tumor:  | 90%   |
| Tumor Hypo/Acellular Stroma:                      | 0%  |
| Tumor Hypercellular Stroma:                       | 10%   |
| Necrosis:   | 0%  |
| Bioanalyzer Ratio (28S/18S):                      | 2.06  |
| CASE Diagnosis (from medical center path report): | Adenocarcinoma of endometrium, endometrioid |
| Age:  | 51  |




[View online »](#)

|                              |  |
|------------------------------|--|
| Gender:                      | Female   |
| Tumor Grade:                 | FIGO G2: Moderately differentiated   |
| Minimum Stage Grouping:      | IB   |
| T (Extent of Primary Tumor): | T1b  |
| N (Lymph node metastasis):   | NX   |
| M (Distant metastasis):      | MX   |
| Storage:                     | Store the total DNA -80°C.   |
| Stability:                   | These products are stable for one year from date of shipping when stored at -80°C without repeated freeze/thaws. |

**Product images:**

4x Image



20x Image